

兔抗 CNN2 多克隆抗体

中文名称：兔抗 CNN2 多克隆抗体

英文名称：Anti-CNN2 rabbit polyclonal antibody

别名：calponin 2

抗原：CNN2

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

相关类别：一抗

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene, which can bind actin, calmodulin, troponin C, and tropomyosin, may function in the structural organization of actin filaments. The encoded protein could play a role in smooth muscle contraction and cell adhesion. Several pseudogenes of this gene have been identified, and are present on chromosomes 1, 2, 3, 6, 9, 11, 13, 15, 16, 21 and 22. Alternative splicing results in multiple transcript variants encoding different isoforms.
Applications:	ELISA, IHC
Name of antibody:	CNN2

Immunogen:	Fusion protein of human CNN2
Full name:	calponin 2
SwissProt:	Q99439
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	20-100

