

## 兔抗 CNPY2 多克隆抗体

中文名称：兔抗 CNPY2 多克隆抗体

英文名称：Anti-CNPY2 rabbit polyclonal antibody

别名：MSAP; TMEM4; ZSIG9; HP10390

抗原：CNPY2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

|                                    |                                                                                                                                                                                       |
|------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Background:</b>                 | Positive regulator of neurite outgrowth by stabilizing myosin regulatory light chain (MRLC). It prevents MIR-mediated MRLC ubiquitination and its subsequent proteasomal degradation. |
| <b>Applications:</b>               | ELISA, WB, IHC                                                                                                                                                                        |
| <b>Name of antibody:</b>           | CNPY2                                                                                                                                                                                 |
| <b>Immunogen:</b>                  | Full length fusion protein                                                                                                                                                            |
| <b>Full name:</b>                  | canopy FGF signaling regulator 2                                                                                                                                                      |
| <b>Synonyms :</b>                  | MSAP; TMEM4; ZSIG9; HP10390                                                                                                                                                           |
| <b>SwissProt:</b>                  | Q9Y2B0                                                                                                                                                                                |
| <b>ELISA Recommended dilution:</b> | 5000-10000                                                                                                                                                                            |

|                                 |                                                  |
|---------------------------------|--------------------------------------------------|
| <b>IHC positive control:</b>    | Human esophagus cancer and human prostate cancer |
| <b>IHC Recommend dilution:</b>  | 25-100                                           |
| <b>WB Predicted band size:</b>  | 21 kDa                                           |
| <b>WB Positive control:</b>     | HEPG2 cell                                       |
| <b>WB Recommended dilution:</b> | 500-2000                                         |



