

兔抗 CNR1 多克隆抗体

中文名称：兔抗 CNR1 多克隆抗体

英文名称：Anti-CNR1 rabbit polyclonal antibody

别名：cannabinoid receptor 1 CB1; CNR; CB-R; CB1A; CB1R; CANN6; CB1K5

相关类别：一抗

储存：冷冻（-20℃）

抗原：CNR1

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes one of two cannabinoid receptors. The cannabinoids, principally delta-9-tetrahydrocannabinol and synthetic analogs, are psychoactive ingredients of marijuana. The cannabinoid receptors are members of the guanine-nucleotide-binding protein (G-protein) coupled receptor family, which inhibit adenylate cyclase activity in a dose-dependent, stereoselective and pertussis toxin-sensitive manner. The two receptors have been found to be involved in the cannabinoid-induced CNS effects (including alterations in mood and cognition) experienced by users of marijuana. Multiple transcript variants encoding tw
--------------------	---

	o different protein isoforms have been described for this gene.
Applications:	ELISA, IHC
Name of antibody:	CNR1
Immunogen:	Synthetic peptide of human CNR1
Full name:	cannabinoid receptor 1
Synonyms:	CB1; CNR; CB-R; CB1A; CB1R; CANN6; CB1K5
SwissProt:	P21554
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer and Human liver cancer
IHC Recommend dilution:	25-100



