

兔抗 CNTN4 多克隆抗体

中文名称：兔抗 CNTN4 多克隆抗体

英文名称：Anti-CNTN4 rabbit polyclonal antibody

别名：contactin 4 AXCAM; BIG-2

抗原：CNTN4

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

相关类别：一抗

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the contactin family of immunoglobulins. Contactins are axon-associated cell adhesion molecules that function in neuronal network formation and plasticity. The encoded protein is a glycosylphosphatidylinositol-anchored neuronal membrane protein that may play a role in the formation of axon connections in the developing nervous system. Deletion or mutation of this gene may play a role in 3p deletion syndrome and autism spectrum disorders. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC

Name of antibody:	CNTN4
Immunogen:	Fusion protein of human CNTN4
Full name:	contactin 4
Synonyms:	AXCAM; BIG-2
SwissProt:	Q8I WV2
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer and Human esophagus cancer
IHC Recommend dilution:	25-100





