

兔抗 FAM186B 多克隆抗体

中文名称：兔抗 FAM186B 多克隆抗体

英文名称：Anti-FAM186B rabbit polyclonal antibody

别名：family with sequence similarity 186 member B; C12orf25

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：FAM186B

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene product is a member of the FAM186 family, however, its exact function is not known. Alternatively spliced transcript variants have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	FAM186B
Immunogen:	Fusion protein of human FAM186B
Full name:	family with sequence similarity 186 member B
Synonyms:	C12orf25
SwissProt:	Q8IYM0
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil

IHC Recommend dilution:	25-100
-------------------------	--------

