



兔抗 FANCD2(Phospho-Ser222) 多克隆抗体

中文名称：兔抗 FANCD2(Phospho-Ser222) 多克隆抗体

英文名称：Anti-FANCD2(Phospho-Ser222) rabbit polyclonal antibody

别 名：FA4; FAD; FACD; FAD2; FA-D2; FANCD

相关类别：一抗

储 存：冷冻 (-20°C) 避光

宿 主：Rabbit

抗 原：FANCD2(Phospho-Ser222)

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The Fanconi anemia complementation group (FANC) currently includes FANCA, FANCB, FANCC, FANCD1 (also called BRC A2), FANCD2, FANCE, FANCF, FANCG, FANCI, FANCJ (also called BRIP1), FANCL, FANCM and FANCN (also called PALB2). The previously defined group FANCH is the same as FANCA. Fanconi anemia is a genetically heterogeneous recessive disorder characterized by cytogenetic instability, hypersensitivity to DNA crosslinking agents, increased chromosomal breakage, and defective DNA repair.
Applications:	WB
Name of antibody:	FANCD2(Phospho-Ser222)

Immunogen:	Peptide sequence around phosphorylation site of serine 222 (G-D-S(p)-Q-H) derived from Human FANCD2.
Full name:	Fanconi anemia, complementation group D2
Synonyms :	FA4; FAD; FACD; FAD2; FA-D2; FANCD
SwissProt:	Q9BXW9
WB Predicted band size:	166 kDa
WB Positive control:	HT29 cells lysates treated with CalyculinA
WB Recommended dilution:	500-1000

