

兔抗 FCGR1A 多克隆抗体

- 中文名称：兔抗 FCGR1A 多克隆抗体
- 英文名称：Anti-FCGR1A rabbit polyclonal antibody
- 别名：Fc fragment of IgG receptor Ia; CD64; FCRI; CD64A; IGFR1
- 相关类别：一抗
- 抗原：FCGR1A
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 反应种属：Human
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a protein that plays an important role in the immune response. This protein is a high-affinity Fc-gamma receptor. The gene is one of three related gene family members located on chromosome 1.
Applications:	ELISA, IHC
Name of antibody:	FCGR1A
Immunogen:	Synthetic peptide of human FCGR1A
Full name:	Fc fragment of IgG receptor Ia
Synonyms:	CD64; FCRI; CD64A; IGFR1
SwissProt:	P12314
ELISA Recommended dilution	5000-10000

n:	
IHC positive control:	Human liver cancer
IHC Recommend dilution:	30-150



