

## CYYR1 抗原（重组蛋白）

中文名称：CYYR1 抗原（重组蛋白）

英文名称：CYYR1 Antigen (Recombinant Protein)

别 名：C21orf95

储 存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 30-154 amino acids of human CYYR1

技术规格：

<b>Full name:</b>	cysteine and tyrosine rich 1
<b>Synonyms:</b>	C21orf95
<b>Swissprot:</b>	Q96J86
<b>Gene Accession:</b>	BC036761
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The function of this protein is unknown.