

兔抗 FEM1A 多克隆抗体

中文名称：兔抗 FEM1A 多克隆抗体

英文名称：Anti-FEM1A rabbit polyclonal antibody

别名：EPRAP

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：FEM1A

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Probable component of an E3 ubiquitin-protein ligase complex, in which it may act as a substrate recognition subunit (By similarity). May participate in antiinflammatory signaling via its interaction with PTGER4.
Applications:	ELISA, IHC
Name of antibody:	FEM1A
Immunogen:	Synthetic peptide of human FEM1A
Full name:	fem-1 homolog A
Synonyms:	EPRAP
SwissProt:	Q9BSK4
ELISA Recommended dilution:	1000-2000

IHC positive control:	Human gasrtic cancer and Human brain
IHC Recommend dilution:	25-100



