

DGCR6L 抗原（重组蛋白）

中文名称：DGCR6L 抗原（重组蛋白）

英文名称：DGCR6L Antigen (Recombinant Protein)

别名：DGCR6

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 21-220 amino acids of human DGCR6L

技术规格：

Full name:	DiGeorge syndrome critical region gene 6 like
Synonyms:	DGCR6
Swissprot:	Q9BY27
Gene Accession:	BC000682
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene, the result of a duplication at this locus, is one of two functional genes encoding nearly identical proteins that have similar expression patterns. The product of this gene is a protein that shares homology with the Drosophila gonadal protein, expressed in gonadal tissues and germ cells, and with the human laminin gamma-1 chain that fu

ctions in cell attachment and migration. This gene is located in a region of chromosome 22 implicated in the DiGeorge syndrome, one facet of a broader collection of anomalies referred to as the CATCH 22 syndrome.