

兔抗 FER(Phospho-Tyr402) 多克隆抗体

中文名称：兔抗 FER(Phospho-Tyr402) 多克隆抗体

英文名称：Anti-FER(Phospho-Tyr402) rabbit polyclonal antibody

别名：FerT; TYK3; Pe1Fe3; Pe1Fe6; PPP1R74; Pe1Fe10; Pe1Fe13

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：FER(Phospho-Tyr402)

反应种属：Human Mouse

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	Fer protein is a member of the FPS/FES family of nontransmembrane receptor tyrosine kinases. It regulates cell-cell adhesion and mediates signaling from the cell surface to the cytoskeleton via growth factor receptors.
Applications:	WB, IHC
Name of antibody:	FER(Phospho-Tyr402)
Immunogen:	Peptide sequence around phosphorylation site of tyrosine 402(V-N-Y(p)-E-E) derived from Human FER .
Full name:	fer (fps/fes related) tyrosine kinase
Synonyms :	FerT; TYK3; Pe1Fe3; Pe1Fe6; PPP1R74; Pe1Fe10; Pe1F

	e13
SwissProt:	P16591
IHC positive control:	Human brain tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	95 kDa
WB Positive control:	Jurkat and COS7 cells lysates
WB Recommended dilution:	500-1000



