

兔抗 FERD3L 多克隆抗体

中文名称：兔抗 FERD3L 多克隆抗体

英文名称：Anti-FERD3L rabbit polyclonal antibody

别名：Fer3 like bHLH transcription factor; PTFB; NATO3; NTWIST; N-TWIST; bHLHa31

相关类别：一抗

储存：冷冻（-20℃）

抗原：FERD3L

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Transcription factor that binds to the E-box and functions as inhibitor of transcription. DNA binding requires dimerization with an E protein. Inhibits transcription activation by ASCL1/MASH1 by sequestering E proteins (By similarity).
Applications:	ELISA, IHC
Name of antibody:	FERD3L
Immunogen:	Synthetic peptide of human FERD3L
Full name:	Fer3 like bHLH transcription factor
Synonyms:	PTFB; NATO3; NTWIST; N-TWIST; bHLHa31
SwissProt:	Q96RJ6

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer and Human cervical cancer
IHC Recommend dilution:	50-300



