

兔抗 FERMT1 多克隆抗体

中文名称：兔抗 FERMT1 多克隆抗体

英文名称：Anti-FERMT1 rabbit polyclonal antibody

别名：fermitin family member 1; URP1; KIND1; DTGCU2; UNC112A; C20orf42

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：FERMT1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the fermitin family, and contains a FERM domain and a pleckstrin homology domain. The encoded protein is involved in integrin signaling and linkage of the actin cytoskeleton to the extracellular matrix. Mutations in this gene have been linked to Kindler syndrome.
Applications:	ELISA, IHC
Name of antibody:	FERMT1
Immunogen:	Synthetic peptide of human FERMT1
Full name:	fermitin family member 1
Synonyms:	URP1; KIND1; DTGCU2; UNC112A; C20orf42

SwissProt:	Q9BQL6
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	25-50

