

兔抗 FES 多克隆抗体

中文名称：兔抗 FES 多克隆抗体

英文名称：Anti-FES rabbit polyclonal antibody

别 名：FPS

相关类别：一抗

抗 原：FES

储 存：冷冻（-20℃）

宿 主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes the human cellular counterpart of a feline sarcoma retrovirus protein with transforming capabilities. The gene product has tyrosine-specific protein kinase activity and that activity is required for maintenance of cellular transformation. Its chromosomal location has linked it to a specific translocation event identified in patients with acute promyelocytic leukemia but it is also involved in normal hematopoiesis as well as growth factor and cytokine receptor signaling. Alternative splicing results in multiple variants encoding different isoforms.
Applications:	ELISA, IHC

Name of antibody:	FES
Immunogen:	Fusion protein of human FES
Full name:	feline sarcoma oncogene
Synonyms :	FPS
SwissProt:	P07332
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human liver cancer and Human esophagus cancer
IHC Recommend dilution:	25-100

