

兔抗 FETUB 多克隆抗体

中文名称：兔抗 FETUB 多克隆抗体

英文名称：Anti-FETUB rabbit polyclonal antibody

别名：16G2; Gugu; IRL685

相关类别：一抗

抗原：FETUB

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the fetuin family, part of the cystatin superfamily of cysteine protease inhibitors. Fetuins have been implicated in several diverse functions, including osteogenesis and bone resorption, regulation of the insulin and hepatocyte growth factor receptors, and response to systemic inflammation. This protein may be secreted by cells.
Applications:	ELISA, WB
Name of antibody:	FETUB
Immunogen:	Fusion protein of human FETUB
Full name:	fetuin B

Synonyms :	16G2; Gugu; IRL685
SwissProt:	Q9UGM5
ELISA Recommended dilution:	2000-5000
WB Predicted band size:	42 kDa
WB Positive control:	231 and A172 cells
WB Recommended dilution:	500-2000

