

## 兔抗 FFAR2 多克隆抗体

中文名称: 兔抗 FFAR2 多克隆抗体

英文名称: Anti-FFAR2 rabbit polyclonal antibody

别名: free fatty acid receptor 2; FFA2R; GPR43

抗原: FFAR2

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the GP40 family of G protein-coupled receptors that are clustered together on chromosome 19. The encoded protein is a receptor for short chain free fatty acids and may be involved in the inflammatory response and in regulating lipid plasma levels.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	FFAR2
<b>Immunogen:</b>	Synthetic peptide of human FFAR2
<b>Full name:</b>	free fatty acid receptor 2
<b>Synonyms:</b>	FFA2R; GPR43

<b>SwissProt:</b>	O15552
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human cervical cancer
<b>IHC Recommend dilution:</b>	25-100

