

兔抗 FGF11 多克隆抗体

中文名称：兔抗 FGF11 多克隆抗体

英文名称：Anti-FGF11 rabbit polyclonal antibody

别名：fibroblast growth factor 11; FHF3; FHF-3; FGF-11

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：FGF11

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the fibroblast growth factor (FGF) family. FGF family members possess broad mitogenic and cell survival activities, and are involved in a variety of biological processes, including embryonic development, cell growth, morphogenesis, tissue repair, tumor growth and invasion. The function of this gene has not yet been determined. The expression pattern of the mouse homolog implies a role in nervous system development. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC

Name of antibody:	FGF11
Immunogen:	Synthetic peptide of human FGF11
Full name:	fibroblast growth factor 11
Synonyms:	FHF3; FHF-3; FGF-11
SwissProt:	Q92914
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human ovarian cancer
IHC Recommend dilution:	25-50

