

兔抗 FLI1 多克隆抗体

中文名称：兔抗 FLI1 多克隆抗体

英文名称：Anti-FLI1 rabbit polyclonal antibody

别名：Fli-1 proto-oncogene, ETS transcription factor; EWSR2; SIC-1; BDPLT21

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：FLI1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a transcription factor containing an ETS DNA-binding domain. The gene can undergo a t(11;22)(q24;q12) translocation with the Ewing sarcoma gene on chromosome 22, which results in a fusion gene that is present in the majority of Ewing sarcoma cases. An acute lymphoblastic leukemia-associated t(4;11)(q21;q23) translocation involving this gene has also been identified. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2012]
Applications:	ELISA, IHC
Name of antibody:	FLI1

Immunogen:	Fusion protein of human FLI1
Full name:	Fli-1 proto-oncogene, ETS transcription factor
Synonyms:	EWSR2; SIC-1; BDPLT21
SwissProt:	Q01543
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	100-200



