

兔抗 FMNL2 多克隆抗体

中文名称：兔抗 FMNL2 多克隆抗体

英文名称：Anti-FMNL2 rabbit polyclonal antibody

别名：formin like 2; FHOD2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：FMNL2

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a formin-related protein. Formin-related proteins have been implicated in morphogenesis, cytokinesis, and cell polarity. Alternatively spliced transcript variants encoding different isoforms have been described but their full-length nature has yet to be determined.
Applications:	ELISA, IHC
Name of antibody:	FMNL2
Immunogen:	Fusion protein of human FMNL2
Full name:	formin like 2
Synonyms:	FHOD2

SwissProt:	Q96PY5
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	50-300

