

HMBOX1 抗原（重组蛋白）

中文名称：HMBOX1 抗原（重组蛋白）

英文名称：HMBOX1 Antigen (Recombinant Protein)

别 名：HOT1; PBHNF; HNF1LA

储 存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 73-272 amino acids of human HMBOX1

技术规格：

Full name:	homeobox containing 1
Synonyms:	HOT1; PBHNF; HNF1LA
Swissprot:	Q6NT76
Gene Accession:	BC009259
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Transcription factor. Ubiquitous. Detected in pancreas, brain, spleen, placenta, prostate, thymus, liver, heart, bone marrow, skeletal muscle, stomach, uterus, testis, kidney, ovary, and colon.