

IFITM2 抗原（重组蛋白）

中文名称：IFITM2 抗原（重组蛋白）

英文名称：IFITM2 Antigen (Recombinant Protein)

别名：1-8D; DSPA2c

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 1-132 amino acids of human IFITM2

技术规格：

Full name:	interferon induced transmembrane protein 2
Synonyms:	1-8D; DSPA2c
Swissprot:	Q01629
Gene Accession:	BC009696
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	IFN-induced antiviral protein which inhibits the entry of viruses to the host cell cytoplasm, permitting endocytosis, but preventing subsequent viral fusion and release of viral contents into the cytosol. Active against multiple viruses, including influenza A virus, SARS coronavirus (SARS-CoV), Marburg virus (MARV), Ebola virus (EBOV), Dengue virus (DENV), West Nile virus (WNV), human immunodeficiency virus type 1 (HIV-1) and vesicular stomatitis virus (VSV). Can inhibit: influenza

	nza virus hemagglutinin protein-mediated viral entry, MARV and EB OV GP1,2-mediated viral entry, SARS-CoV S protein-mediated viral e ntry and VSV G protein-mediated viral entry. Induces cell cycle arres t and mediates apoptosis by caspase activation and in p53-indepen dent manner.
--	---