

兔抗 FNTA 多克隆抗体

中文名称：兔抗 FNTA 多克隆抗体

英文名称：Anti-FNTA rabbit polyclonal antibody

别名：FPTA; PTAR2; PGGT1A

相关类别：一抗

抗原：FNTA

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Prenyltransferases can attach either a farnesyl group or a geranylgeranyl group in thioether linkage to the cysteine residue of proteins with a C-terminal CAAX box. CAAX geranylgeranyltransferase and CAAX farnesyltransferase are heterodimers that share the same alpha subunit but have different beta subunits. This gene encodes the alpha subunit of these transferases. Alternative splicing results in multiple transcript variants. Related pseudogenes have been identified on chromosomes 11 and 13.
Applications:	ELISA, IHC
Name of antibody:	FNTA

Immunogen:	Full length fusion protein
Full name:	farnesyltransferase, CAAX box, alpha
Synonyms:	FPTA; PTAR2; PGGT1A
SwissProt:	P49354
IHC positive control:	Human thyroid cancer and human liver cancer
IHC Recommend dilution:	25-100



