

LDHAL6B 抗原（重组蛋白）

中文名称：LDHAL6B 抗原（重组蛋白）

英文名称：LDHAL6B Antigen (Recombinant Protein)

别 名：LDHL; LDH6B; LDHAL6

储 存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 182-381 amino acids of human LDHAL6B

技术规格：

Full name:	lactate dehydrogenase A like 6B
Synonyms:	LDHL; LDH6B; LDHAL6
Swissprot:	Q9BYZ2
Gene Accession:	BC022034
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	LDHAL6B belongs to the LDH/MDH superfamily. LDH family. The specific function of the protein remains unknown.