

兔抗 FOXH1 多克隆抗体

中文名称：兔抗 FOXH1 多克隆抗体

英文名称：Anti-FOXH1 rabbit polyclonal antibody

别名：forkhead box H1; FAST1; FAST-1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：FOXH1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	FOXH1 encodes a human homolog of Xenopus forkhead activin signal transducer-1. FOXH1 protein binds SMAD2 and activates an activin response element via binding the DNA motif TGT(G/T)(T/G)ATT.
Applications:	ELISA, IHC
Name of antibody:	FOXH1
Immunogen:	Synthetic peptide of human FOXH1
Full name:	forkhead box H1
Synonyms:	FAST1; FAST-1
SwissProt:	O75593
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer

IHC Recommend dilution:

50-100

