

兔抗 FOXI1 多克隆抗体

中文名称：兔抗 FOXI1 多克隆抗体

英文名称：Anti-FOXI1 rabbit polyclonal antibody

别名：forkhead box I1; HFH3; FKH10; HFH-3; FKHL10; FREAC6; FREAC-6

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：FOXI1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene belongs to the forkhead family of transcription factors, which is characterized by a distinct forkhead domain. This gene may play an important role in the development of the cochlea and vestibulum, as well as in embryogenesis. The encoded protein has been found to be required for the transcription of four subunits of a proton pump found in the inner ear, the kidney, and the epididymis. Mutations in this gene have been associated with deafness, autosomal recessive 4.
Applications:	ELISA, WB, IHC
Name of antibody:	FOXI1

Immunogen:	Synthetic peptide of human FOXI1
Full name:	forkhead box I1
Synonyms:	HFH3; FKH10; HFH-3; FKHL10; FREAC6; FREAC-6
SwissProt:	Q12951
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human ovarian cancer and Human tonsil
IHC Recommend dilution:	50-200
WB Predicted band size:	41 kDa
WB Positive control:	Human heart tissue lysate
WB Recommended dilution:	500-2000



