

## 兔抗 FOXO1 (Phospho-Ser256)多克隆抗体

中文名称：兔抗 FOXO1 (Phospho-Ser256)多克隆抗体

英文名称：Anti-FOXO1 (Phospho-Ser256) rabbit polyclonal antibody

别名：FKH1; FKHR; FOXO1A

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene belongs to the forkhead family of transcription factors which are characterized by a distinct forkhead domain. The specific function of this gene has not yet been determined; however, it may play a role in myogenic growth and differentiation. Translocation of this gene with PAX3 has been associated with alveolar rhabdomyosarcoma.
<b>Applications:</b>	WB, IHC, IF
<b>Name of antibody:</b>	FOXO1 (Phospho-Ser256)
<b>Immunogen:</b>	Synthetic peptide of human FOXO1 (Phospho-Ser256)
<b>Full name:</b>	forkhead box O1 (Phospho-Ser256)
<b>Synonyms :</b>	FKH1; FKHR; FOXO1A
<b>SwissProt:</b>	Q12778
<b>IHC positive control:</b>	Human breast carcinoma

<b>IHC Recommend dilution:</b>	50-100
<b>WB Predicted band size:</b>	70 kDa
<b>WB Positive control:</b>	293 cells untreated or treated with serum
<b>WB Recommended dilution:</b>	500-1000
<b>IF Positive control:</b>	MCF7 cells
<b>IF Recommended dilution</b>	100-200



