

MPV17 抗原（重组蛋白）

中文名称：MPV17 抗原（重组蛋白）

英文名称：MPV17 Antigen (Recombinant Protein)

别名：SYM1; MTDPS6

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 39-176 amino acids of human MPV17

技术规格：

Full name:	MPV17, mitochondrial inner membrane protein
Synonyms:	SYM1; MTDPS6
Swissprot:	P39210
Gene Accession:	BC001115
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a mitochondrial inner membrane protein that is implicated in the metabolism of reactive oxygen species. Mutations in this gene have been associated with the hepatocerebral form of mitochondrial DNA depletion syndrome (MDDS).