

## MRPL50 抗原(重组蛋白)

中文名称: MRPL50 抗原(重组蛋白)

- 英文名称: MRPL50 Antigen (Recombinant Protein)
- 别名: MRP-L50
- 储存: 冷冻(-20℃)
- 相关类别: 抗原

概述:

Fusion protein corresponding to a region derived from 1-158 amino acids of human MRPL50

## 技术规格:

Full name:	mitochondrial ribosomal protein L50
Synonyms:	MRP-L50
Swissprot:	Q8N5N7
Gene Accession:	BC032008
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Mammalian mitochondrial ribosomal proteins are encoded by nucle ar genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S su bunit and a large 39S subunit. They have an estimated 75% protei n to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitor ibosomes and prokaryotic ribosomes is that the latter contain a 5S



rRNA. Among different species, the proteins comprising the mitorib
osome differ greatly in sequence, and sometimes in biochemical pr
operties, which prevents easy recognition by sequence homology. T
his gene encodes a putative 39S subunit protein and belongs to t
he L47P ribosomal protein family. Pseudogenes corresponding to t
his gene are found on chromosomes 2p, 2q, 5p, and 10q.