

NADK 抗原（重组蛋白）

中文名称：NADK 抗原（重组蛋白）

英文名称：NADK Antigen (Recombinant Protein)

别 名：dJ283E3.1

储 存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 247-446 amino acids of human NADK

技术规格：

Full name:	NAD kinase
Synonyms:	dJ283E3.1
Swissprot:	O95544
Gene Accession:	BC001709
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	NADK catalyzes the transfer of a phosphate group from ATP to NAD to generate NADP, which in its reduced form acts as an electron donor for biosynthetic reactions (Lerner et al., 2001 [PubMed 11594753]).