

兔抗 FOXQ1 多克隆抗体

中文名称：兔抗 FOXQ1 多克隆抗体

英文名称：Anti-FOXQ1 rabbit polyclonal antibody

别名：forkhead box Q1; HFH1

抗原：FOXQ1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	FOXQ1 is a member of the FOX gene family, which is characterized by a conserved 110-amino acid DNA-binding motif called the forkhead or winged helix domain. FOX genes are involved in embryonic development, cell cycle regulation, tissue-specific gene expression, cell signaling, and tumorigenesis (Bieller et al., 2001 [PubMed 11747606]).
Applications:	ELISA, WB, IHC
Name of antibody:	FOXQ1
Immunogen:	Synthetic peptide of human FOXQ1
Full name:	forkhead box Q1

Synonyms:	HFH1
SwissProt:	Q9C009
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human gastric cancer and Human ovarian cancer
IHC Recommend dilution:	25-50
WB Predicted band size:	42 kDa
WB Positive control:	Mouse kidney tissue lysate
WB Recommended dilution:	500-2000



