

兔抗 FPR2 多克隆抗体

中文名称：兔抗 FPR2 多克隆抗体

英文名称：Anti-FPR2 rabbit polyclonal antibody

别名：ALXR, HM63, FMLPX, FPR2A, FPRH1, FPRH2, FPRL1, LXA4R, FMLP-R-II

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：FPR2

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Low affinity receptor for N-formyl-methionyl peptides, which are powerful neutrophils chemotactic factors. Binding of FMLP to the receptor causes activation of neutrophils. This response is mediated via a G-protein that activates a phosphatidylinositol-calcium second messenger system. The activation of LXA4R could result in an anti-inflammatory outcome counteracting the actions of proinflammatory signals such as LTB4
Applications:	ELISA, IHC
Name of antibody:	FPR2

Immunogen:	Synthetic peptide of human FPR2
Full name:	formyl peptide receptor 2
Synonyms :	ALXR, HM63, FMLPX, FPR2A, FPRH1, FPRH2, FPRL1, LXA4R, FMLP-R-II
SwissProt:	P25090
ELISA Recommended dilution:	1000-10000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	20-100

