

兔抗 FRS3 多克隆抗体

中文名称：兔抗 FRS3 多克隆抗体

英文名称：Anti-FRS3 rabbit polyclonal antibody

别名：SNT2; FRS2B; SNT-2; FRS2beta; FRS2-beta

储存：冷冻（-20℃）避光

抗原：FRS3

宿主：Rabbit

反应种属：Human

相关类别：一抗

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	his gene encodes a substrate for the fibroblast growth factor receptor. The encoded protein is found in the peripheral plasma membrane and links fibroblast growth factor receptor stimulation to activators of Ras. The encoded protein down-regulates extracellular regulated kinase 2 through direct binding.
Applications:	WB
Name of antibody:	FRS3
Immunogen:	Synthesized peptide derived from internal of human FRS3.
Full name:	fibroblast growth factor receptor substrate 3

Synonyms :	SNT2; FRS2B; SNT-2; FRS2beta; FRS2-beta
SwissProt:	O43559
WB Predicted band size:	54 kDa
WB Positive control:	COS-7 cells lysate
WB Recommended dilution:	500-3000

