

兔抗 FSD1 多克隆抗体

中文名称：兔抗 FSD1 多克隆抗体

英文名称：Anti-FSD1 rabbit polyclonal antibody

别名：MIR1; GLFND

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：FSD1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a centrosome associated protein that is characterized by an N-terminal coiled-coil region downstream of B-box (BBC) domain, a central fibronectin type II I domain, and a C-terminal repeats in splA and RyR (SPRY) domain. The encoded protein associates with a subset of microtubules and may be involved in the stability and organization of microtubules during cytokinesis.
Applications:	ELISA, WB, IHC
Name of antibody:	FSD1
Immunogen:	Fusion protein of human FSD1
Full name:	fibronectin type III and SPRY domain containing 1

Synonyms:	MIR1; GLFND
SwissProt:	Q9BTV5
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human cervical cancer and Human colorectal cancer
IHC Recommend dilution:	20-100
WB Predicted band size:	56 kDa
WB Positive control:	SP20 cell lysate
WB Recommended dilution:	200-1000



