

兔抗 FSTL5 多克隆抗体

中文名称：兔抗 FSTL5 多克隆抗体

英文名称：Anti-FSTL5 rabbit polyclonal antibody

别名：follistatin like 5

抗原：FSTL5

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Follistatin-like 5 (FSTL5), a member of the follistatin family of genes, encodes a secretory glycoprotein. Previous study revealed that it might play a suppressive role in hepatocellular carcinoma (HCC). However, its clinical significances, biological functions and molecular mechanisms in HCC development are poorly understood.
Applications:	ELISA, IHC
Name of antibody:	FSTL5
Immunogen:	Synthetic peptide of human FSTL5
Full name:	follistatin like 5
SwissProt:	Q8N475

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	50-300

