

兔抗 FTSJ3 多克隆抗体

中文名称：兔抗 FTSJ3 多克隆抗体

英文名称： Anti-FTSJ3 rabbit polyclonal antibody

别名： FtsJ homolog 3; SPB1; EPCS3

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： FTSJ3

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Although the function of this gene is not known, the existence of this gene is supported by mRNA and EST data. A possible function of the encoded protein can be inferred from amino acid sequence similarity to the E.coli FtsJ protein and to a mouse protein possibly involved in embryogenesis.
Applications:	ELISA, IHC
Name of antibody:	FTSJ3
Immunogen:	Fusion protein of human FTSJ3
Full name:	FtsJ homolog 3
Synonyms:	SPB1; EPCS3

SwissProt:	Q8IY81
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human ovarian cancer
IHC Recommend dilution:	30-150

