

兔抗 GLYAT 多克隆抗体

中文名称：兔抗 GLYAT 多克隆抗体

英文名称：Anti-GLYAT rabbit polyclonal antibody

别名：glycine-N-acyltransferase; CAT; GAT; ACGNAT

相关类别：一抗

储存：冷冻（-20℃）

抗原：GLYAT

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The glycine-N-acyltransferase protein conjugates glycine with acyl-CoA substrates in the mitochondria. The protein is thought to be important in the detoxification of endogenous and xenobiotic acyl-CoA's. Two transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, WB, IHC
Name of antibody:	GLYAT
Immunogen:	Fusion protein of human GLYAT
Full name:	glycine-N-acyltransferase
Synonyms:	CAT; GAT; ACGNAT

SwissProt:	Q6IB77
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human ovarian cancer
IHC Recommend dilution:	50-300
WB Predicted band size:	34 kDa
WB Positive control:	Mouse kidney tissue and Human kidney tissue lysates
WB Recommended dilution:	500-2000



