

兔抗 GNAO1 多克隆抗体

中文名称：兔抗 GNAO1 多克隆抗体

英文名称：Anti-GNAO1 rabbit polyclonal antibody

别名：GNAO; EIEE17; G-ALPHA-o

相关类别：一抗

克隆类型：rabbit polyclonal

抗原：GNAO1

储存：冷冻（-20℃）避光

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

技术规格

Background:	Guanine nucleotide-binding proteins (G proteins) are involved as modulators or transducers in various transmembrane signaling systems. The G(o) protein function is not clear. Stimulated by RGS14.
Applications:	WB, IHC
Name of antibody:	GNAO1
Immunogen:	Fusion protein of human GNAO1
Full name:	guanine nucleotide binding protein (G protein), alpha activating activity polypeptide O
Synonyms :	GNAO; EIEE17; G-ALPHA-o
SwissProt:	P09471

IHC positive control:	Human brain
IHC Recommend dilution:	25-100
WB Predicted band size:	40 kDa
WB Positive control:	Mouse brain tissue, hela cell and human fetal brain tissue
WB Recommended dilution:	500-2000



