

兔抗 GPAM 多克隆抗体

中文名称：兔抗 GPAM 多克隆抗体

英文名称：Anti-GPAM rabbit polyclonal antibody

别名：GPAT; GPAT1

相关类别：一抗

抗原：GPAM

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a mitochondrial enzyme which prefers saturated fatty acids as its substrate for the synthesis of glycerolipids. This metabolic pathway's first step is catalyzed by the encoded enzyme. Two forms for this enzyme exist, one in the mitochondria and one in the endoplasmic reticulum. Two alternatively spliced transcript variants have been described for this gene.
Applications:	ELISA, IHC
Name of antibody:	GPAM
Immunogen:	Fusion protein of human GPAM
Full name:	glycerol-3-phosphate acyltransferase, mitochondrial

Synonyms :	GPAT; GPAT1
SwissProt:	Q9HCL2
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human ovarian cancer and human gastric cancer
IHC Recommend dilution:	100-300

