

## 兔抗 GPANK1 多克隆抗体

中文名称: 兔抗 GPANK1 多克隆抗体

英文名称: Anti-GPANK1 rabbit polyclonal antibody

别名: G5; BAT4; D6S54E; ANKRD59; GPATCH10

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: GPANK1

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

相关类别: 一抗

### 技术规格

<b>Background:</b>	This gene is located in a cluster of HLA-B-associated transcripts, which is included in the human major histocompatibility complex III region. This gene encodes a protein which is thought to play a role in immunity. Multiple alternatively spliced variants, encoding the same protein, have been identified.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	GPANK1
<b>Immunogen:</b>	Fusion protein of human GPANK1
<b>Full name:</b>	G-patch domain and ankyrin repeats 1
<b>Synonyms:</b>	G5; BAT4; D6S54E; ANKRD59; GPATCH10

<b>SwissProt:</b>	O95872
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	25-100



