

兔抗 GPATCH1 多克隆抗体

中文名称：兔抗 GPATCH1 多克隆抗体

英文名称：Anti-GPATCH1 rabbit polyclonal antibody

别名：G-patch domain containing 1; ECGP; GPATC1

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：GPATCH1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The function of this protein remains unknown.
Applications:	ELISA, WB, IHC
Name of antibody:	GPATCH1
Immunogen:	Synthetic peptide of human GPATCH1
Full name:	G-patch domain containing 1
Synonyms:	ECGP; GPATC1
SwissProt:	Q9BRR8
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and Human esophagus cancer
IHC Recommend dilution:	20-100
WB Predicted band size:	103 kDa

WB Positive control:	A549 cell lysate
WB Recommended dilution:	200-1000



