

TTC16 抗原（重组蛋白）

中文名称：TTC16 抗原（重组蛋白）

英文名称：TTC16 Antigen (Recombinant Protein)

别名：tetratricopeptide repeat domain 16

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 674-873 amino acids of human TTC16

技术规格：

Full name:	tetratricopeptide repeat domain 16
Swissprot:	Q8NEE8
Gene Accession:	BC031281
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The function of this protein remains unknown.