

兔抗 GPR21 多克隆抗体

中文名称: 兔抗 GPR21 多克隆抗体

英文名称: Anti-GPR21 rabbit polyclonal antibody

别名: G protein-coupled receptor 21

抗原: GPR21

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a member of the G-protein-coupled receptor 1 family. G-protein coupled receptors are membrane proteins which activate signaling cascades as a response to extracellular stress. The encoded protein activates a Gq signal transduction pathway which mobilizes calcium.
Applications:	ELISA, IHC
Name of antibody:	GPR21
Immunogen:	Synthetic peptide of human GPR21
Full name:	G protein-coupled receptor 21
SwissProt:	Q99679

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer and Human cervical cancer
IHC Recommend dilution:	30-150



