

兔抗 HBP1 多克隆抗体

- 中文名称：兔抗 HBP1 多克隆抗体
- 英文名称：Anti-HBP1 rabbit polyclonal antibody
- 相关类别：一抗
- 储 存：冷冻（-20℃）避光
- 抗 原：HBP1
- 宿 主：Rabbit
- 反应种属：Human Mouse Rat
- 标记物：Unconjugate
- 克隆类型：Unconjugate

技术规格

Background:	Transcriptional repressor that binds to the promoter region of target genes. Plays a role in the regulation of the cell cycle and of the Wnt pathway. Binds preferentially to the sequence 5'-TTCATTCATTCA-3'. Binding to the H1F0 promoter is enhanced by interaction with RB1. Disrupts the interaction between
--------------------	---

	en DNA and TCF4.
Applications:	WB
Name of antibody:	HBP1
Immunogen:	Synthesized peptide derived from Internal of human HBP1.
Full name:	HMG-box transcription factor 1
Synonyms :	0
SwissProt:	O60381
WB Predicted band size:	58 kDa
WB Positive control:	COS-7 cells lysate
WB Recommended dilution:	500-3000

