

兔抗 GPR27 多克隆抗体

中文名称: 兔抗 GPR27 多克隆抗体

英文名称: Anti-GPR27 rabbit polyclonal antibody

别名: G protein-coupled receptor 27; SREB1

储存: 冷冻 (-20℃)

抗原: GPR27

宿主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	GPR27 is a member of the G protein-coupled receptors (GPCRs), a large family of receptors that have a similar structure characterized by 7 transmembrane domains. Activation of GPCRs by extracellular stimuli such as neurotransmitters, hormones, or light induces an intracellular signaling cascade mediated by heterotrimeric GTP-binding proteins, or G proteins.
Applications:	ELISA, IHC
Name of antibody:	GPR27
Immunogen:	Synthetic peptide of human GPR27
Full name:	G protein-coupled receptor 27

Synonyms:	SREB1
SwissProt:	Q9NS67
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer and Human tonsil
IHC Recommend dilution:	30-150



