

## 兔抗 MAGEB5 多克隆抗体

中文名称：兔抗 MAGEB5 多克隆抗体

英文名称：Anti-MAGEB5 rabbit polyclonal antibody

别名：CT3.3; MAGE-B5

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：MAGEB5

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

#### Background:

The melanoma-associated antigen (MAGE) family consists of a number of antigens recognized by cytotoxic T lymphocytes. The MAGE genes were initially isolated from different kinds of tumors and, based on their virtually exclusive tumor-specific expression in adult tissues, they have been used as targets for cancer immunotherapy. MAGE genes encode for tumor-rejection antigens that are expressed in tumors of different histologic types and in normal testis and placenta. MAGE-B5 (melanoma-associated antigen B5), also known as CT3.3 (cancer/testis antigen 3.3), is a 275 amino acid protein that con

	<p>tains one MAGE domain and is expressed in testis. The gene encoding MAGE-B5 maps to human chromosome X, which consists of about 153 million base pairs and nearly 1,000 genes. Color blindness, hemophilia and Du chenne muscular dystrophy are well known X chromosome-linked conditions which affect males more frequently, as males carry a single X chromosome.</p>
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	MAGEB5
<b>Immunogen:</b>	Synthetic peptide of human MAGEB5
<b>Full name:</b>	MAGE family member B5
<b>Synonyms:</b>	CT3.3; MAGE-B5
<b>SwissProt:</b>	Q9BZ81
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human thyroid cancer and human liver cancer
<b>IHC Recommend dilution:</b>	20-100



