

兔抗 GPR142 多克隆抗体

中文名称: 兔抗 GPR142 多克隆抗体

英文名称: Anti-GPR142 rabbit polyclonal antibody

别名: G protein-coupled receptor 142; PGR2; GPRg1b

储存: 冷冻 (-20℃)

抗原: GPR142

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	GPR142 is a member of the rhodopsin family of G protein-coupled receptors (GPRs).
Applications:	ELISA, WB
Name of antibody:	GPR142
Immunogen:	Synthetic peptide of human GPR142
Full name:	G protein-coupled receptor 142
Synonyms:	PGR2; GPRg1b
SwissProt:	Q7Z601
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	51 kDa
WB Positive control:	Hela cell lysate

WB Recommended dilution:

500-2000





