

TUBD1 抗原（重组蛋白）

中文名称：TUBD1 抗原（重组蛋白）

英文名称： TUBD1 Antigen (Recombinant Protein)

别 名： TUBD

储 存： 冷冻（-20℃）

相关类别： 抗原

概述：

Fusion protein corresponding to a region derived from 1-200 amino acids of human TUBD1

技术规格：

Full name:	tubulin delta 1
Synonyms:	TUBD
Swissprot:	Q9UJT1
Gene Accession:	BC000258
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	In the elongating spermatid it is associated with the manchette, a specialized microtubule system present during reshaping of the sperm head.