

UBAC1 抗原（重组蛋白）

中文名称：UBAC1 抗原（重组蛋白）

英文名称：UBAC1 Antigen (Recombinant Protein)

别名：KPC2; GBDR1; UBADC1

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 206-405 amino acids of human UBAC1

技术规格：

Full name:	UBA domain containing 1
Synonyms:	KPC2; GBDR1; UBADC1
Swissprot:	Q9BSL1
Gene Accession:	BC004967
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Non-catalytic subunit of the KPC complex that acts as E3 ubiquitin-protein ligase. Required for poly-ubiquitination and proteasome-mediated degradation of CDKN1B during G1 phase of the cell cycle.