

## UNKL 抗原（重组蛋白）

中文名称：UNKL 抗原（重组蛋白）

英文名称：UNKL Antigen (Recombinant Protein)

别名：ZC3H5L; C16orf28; ZC3HDC5L

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 1-277 amino acids of human UNKL

技术规格：

<b>Full name:</b>	unkempt family like zinc finger
<b>Synonyms:</b>	ZC3H5L; C16orf28; ZC3HDC5L
<b>Swissprot:</b>	Q9H9P5
<b>Gene Accession:</b>	BC011924
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a RING finger protein that may function in Rac signaling. It can bind to Brg/Brm-associated factor 60b and can promote its ubiquitination. Alternative splicing of this gene results in multiple transcript variants.