

兔抗 GPR151 多克隆抗体

中文名称: 兔抗 GPR151 多克隆抗体

英文名称: Anti-GPR151 rabbit polyclonal antibody

别名: G protein-coupled receptor 151; GPCR; PGR7; GALR4; GALRL; GPCR-2037

储存: 冷冻 (-20℃)

抗原: GPR151

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is similar to the galanin receptor subfamily of G protein-coupled receptors . The encoded protein is found predominantly in the central nervous system.
Applications:	ELISA, IHC
Name of antibody:	GPR151
Immunogen:	Synthetic peptide of human GPR151
Full name:	G protein-coupled receptor 151
Synonyms:	GPCR; PGR7; GALR4; GALRL; GPCR-2037
SwissProt:	Q8TDV0
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human cervical cancer and Human colorectal cancer
IHC Recommend dilution:	40-200



